

Product of the Month

SERVA
■ serving scientists ■



October

Get **20%** discount* now!

*For all SERVA customers within the EU. Outside Germany and Austria online shop orders only. Can not be combined with other discounts, promotions etc.

Protease Inhibitor Mixes

Effective protection against protein degradation

Optimized for protein isolation from various sources

During preparation of cell extracts, proteases are inevitably released from bacteria, yeast, tissue or cell culture.

SERVA Protease Inhibitor Mixes provide optimal protection of proteins against proteolytic degradation.

Features

- Easy handling
- High solubility
- Application-optimized
- One vial yields 1 ml of a 100-fold concentrate



Protease Inhibitor Mix G

- For general applications
- Contains 5 protease inhibitors
- Soluble in water

Protease Inhibitor Mix M

- For extracts from mammalian cells
- Contains 6 protease inhibitors
- DMSO for resuspension included

Protease Inhibitor Mix FY

- For fungi and yeast extracts
- Contains 4 protease inhibitors
- DMSO for resuspension included

Protease Inhibitor Mix HP

- Protection of recombinant polyHis-tagged proteins expressed in *E. coli*
- Contains 4 protease inhibitors
- Soluble in water

Ordering information

Product	Cat. no.	Size	List price	Special price*
Protease Inhibitor Mix G	39101.02	5 vials	127.00 €	101.60 €
	39101.03	10 vials	224.00 €	179.20 €
Protease Inhibitor Mix M	39102.02	5 vials	207.00 €	165.60 €
	39102.03	10 vials	351.00 €	280.80 €
Protease Inhibitor Mix FY	39104.02	5 vials	172.00 €	137.60 €
	39104.03	10 vials	314.00 €	251.20 €
Protease Inhibitor Mix HP	39106.02	5 vials	127.00 €	101.60 €
	39106.03	10 vials	224.00 €	179.20 €

* The special prices are valid until October 31st, 2021.

For further information, please contact:

Dr. Judith Koch
Tel.: +49 6221 13840-44
E-Mail: Tech.Service@serva.de

SERVA Electrophoresis GmbH
Carl-Benz-Str. 7
D-69115 Heidelberg • Germany
www.serva.de