


# Horizontal Protein Electrophoresis



## Excellent Gel Kits

- Significantly higher band sharpness compared to vertical SDS PAGE
- Behaves in separation performance and handling like ExcelGel® and can be used as a direct substitute
- Available with non-fluorescent (NF) or standard support foil depending on the staining method
- Format 260 mm x 125 mm x 0.43 mm

**- 25%**

Cat. no.	Product	Size	Price	Promotion
43421.01	Excellent Gel Kit 12.5 % for 1D SDS PAGE, contains 4 gels and buffer kit	1 kit	222,00 €	166,50 €
43422.01	Excellent Gel Kit 7.5 % for 1D SDS PAGE, contains 4 gels and buffer kit	1 kit	222,00 €	166,50 €
43423.01	Excellent Gel Kit NF 12.5 % for 1D SDS PAGE, contains 4 gels and buffer kit	1 kit	222,00 €	166,50 €
43424.01	Excellent Gel Kit NF 7.5 % for 1D SDS PAGE, contains 4 gels and buffer kit	1 kit	222,00 €	166,50 €
43425.01	Excellent Gel Kit NF 15 % for 1D SDS PAGE, contains 4 gels and buffer kit	1 kit	222,00 €	166,50 €
43426.01	Excellent Gel Kit NF 15 % for 1D SDS PAGE, contains 4 gels and buffer kit	1 kit	222,00 €	166,50 €


## Ready-to-use Gels for the PhastSystem™

- Optionally for IEF or SDS PAGE
- Manufactured in Heidelberg

**New!**

**- 25%**

Cat. no.	Product	Content	Size	Price	Promotion
43366.01	SERVA IEF Gel 3 - 9 for PhastSystem™ Size 50 mm x 42 mm x 0.43 mm	10 gels	10 gels	175,00 €	131,25 €
43503.01	SERVA SDS Gel 8 - 25 % Kit for PhastSystem™ Size 50 mm x 42 mm x 0.43 mm	10 gels, electrode wicks, cathode and anode buffer	1 kit	225,00 €	168,75 €



## 2D HPE™ BlotGel NF 12.5 %

- For Western blots after high-resolution, horizontal 2D electrophoresis
- Optimised gels for semi-dry blotting due to non-covalent adhesion to carrier film
- Size 250 mm x 110 mm x 0.65 mm

**New!**

**- 25%**

Cat. no.	Product	Content	Size	Price	Promotion
43429.01	2D HPE™ Triple BlotGel NF 12.5 % Kit	4 gels for three 7 cm IPG strips each, buffer, electrode wicks	1 kit	330,00 €	247,50 €
43430.01	2D HPE™ Double BlotGel NF 12.5 % Kit	4 gels for two 11 cm IPG strips each, buffer, electrode wicks	1 kit	330,00 €	247,50 €
43432.01	2D HPE™ Large BlotGel NF 12.5 % Kit	4 gels for one 24 cm IPG strip each, buffer, electrode wicks	1 kit	400,00 €	300,00 €




## Staining kits

- Optimized protocols for staining thin, horizontal gels on carrier film

**- 25%**

Cat. no.	Product	Size	Price	Promotion
43395.01	SERVA HPE™ Silver Staining Kit	1 kit	122,00 €	91,50 €
43396.01	SERVA HPE™ Coomassie® Staining Kit	1 kit	96,00 €	72,00 €



## HPE™ BlueHorizon™

- Optimized performance with all types of flatbed gel electrophoresis
- Highest resolution due to premium quality ceramic cooling plate
- Workshops and IQ/OQ/PQ on request

Cat. no.	Product	Size	Price	Promotion
HPE-BH	HPE™ BlueHorizon™	1 piece	5.700,00 €	on request



The power supply unit shown is not included in the scope of delivery of the HPE™ BlueHorizon™ (HPE-BH).

EXCELGel™ and PhastSystem™ = Trademarks from Cytiva/former GE Healthcare

# SERVA

Promotion is valid until June 30<sup>th</sup>, 2021  
 How to order: phone: +49 6221 13840-47, fax: -10, e-mail: order@serva.de or visit the SERVA online shop at www.serva.de.  
 For SERVA General Conditions of Sale and Supply (GCSS) please refer to www.serva.de/gcss.