


Western Blotting

SERVAGe™ Neutral HSE

- HSE - High Speed Electrophoresis: Gel run completed in 20 minutes at 400 V
- Separation range from 6.5 kDa to 200 kDa
- For standard SDS Tris/Glycine buffer system according to Laemmli
- 15 months shelf life, in cassettes 10 cm x 10 cm x 0.7 cm

- 25%

Cat. no.	Product	Wells	Size	Price	Promotion
43246.01	SERVAGe™ Neutral HSE gel	10	10 Gele	132,00 €	99,00 €
43245.01		12		132,00 €	99,00 €
43249.01		15		132,00 €	99,00 €
43247.01		2D		132,00 €	99,00 €



Xpress Blotting Kit

- Optimized for fast and efficient semi-dry transfer in only 15 min
- Blotting fleece for particularly even transfer

- 30%

Cat. no.	Product	Content	Size	Price	Promotion
42661.01	Xpress Blotting Buffer (10x)	Xpress Blotting Buffer (10x)	1 l	76,00 €	53,20 €
42662.01	Xpress Blotting Kit	Xpress Blotting Buffer (10x), Blotting Fleece (80 mm x 85 mm)	1 kit (for 10 gels)	61,00 €	42,70 €
42663.01	Xpress NC Blotting Kit	Xpress Blotting Buffer (10x), Blotting Fleece, NC Sheets (80 x 85 mm, 0.2 µm)	1 kit (for 10 gels)	139,00 €	97,30 €
42664.01	Xpress PVDF Blotting Kit	Xpress Blotting Puffer (10x), Blotting Fleece, PVDF Sheets (80 x 85 mm, 0.2 µm)	1 kit (for 10 gels)	190,00 €	133,00 €



BlueBlock PF

- Protein-free blocking solution
- Low background with amplified signal intensity

- 30%

Cat. no.	Product	Size	Price	Promotion
42591.01	BlueBlock PF (10x)	250 ml	44,00 €	30,80 €
42591.02		1 l	136,00 €	95,20 €



BlueClear SB Antibody Stripping Buffer

- Efficiently removes even highly specific antibodies after CL detection
- Without β-ME/DTT, suitable for NC and PVDF membranes

- 25%

Cat. no.	Product	Size	Price	Promotion
42599.01	BlueClear SB Antibody Stripping Buffer	250 ml	40,00 €	30,00 €
42599.02		1 l	99,00 €	74,25 €




BlueBlot Semi Dry Blotter

- Platinum covered steel net cathode and anode electrodes for homogeneous protein transfer
- Spring mounted anode for flexible applications
- Stable acryl housing, device and electrodes are very easy to clean

- 20%

Cat. no.	Product	Blotting area	Size	Price	Promotion
BB-SD11	BlueBlot Semi Dry Blotter SD11	11 cm x 11 cm	1 piece	2.025,00 €	1.620,00 €
BB-SD17	BlueBlot Semi Dry Blotter SD17	17 cm x 17 cm	1 piece	2.500,00 €	2.000,00 €
BB-SD26	BlueBlot Semi Dry Blotter SD26	24 cm x 26 cm	1 piece	3.230,00 €	2.584,00 €



SERVA

Promotion is valid until June 30th, 2021
 How to order: phone: +49 6221 13840-47, fax: -10, e-mail: order@serva.de or visit the SERVA online shop at www.serva.de.
 For SERVA General Conditions of Sale and Supply (GCSS) please refer to www.serva.de/gcss.