

Anwendungen (1)	Membranproteine/ Proteomforschung	Rekonstitution von Membran- proteinen	Rezeptor Forschung	Enzymologie	Photosynthese- Forschung	Reverse Mizellen- löslichkeit	Lipid Raft Forschung	DNA- Isolierung	RNA- isolierung
Brij 35	X			X					
CHAPS	X		X	X	X				
CTAP						X		X	X
Na-Cholsäure		X	X	X					
Na- Deoxycholsäure			X						
Digitonin	X		X	X	X				
Na-Docusat						X			
DM	X		X	X	X				
SDS	X			X					
L-Sarcosin	X			X				X	X
Lubrol	X			X			X		
OG	X	X	X	X	X				
Saponin									
SB 10	X			X					
SB 12					X				
Synperonic F68									
Synperonic F108									
Triton® X-100	X		X	X	X		X		
Tween 20				X					
Tween 80				X					
ASB-14/16	X								

Anwendungen (2)	Enzym- Assays	Immuno- Assays	Chemilumineszenz- Analyse/Fluorometrie	Chromatographie	Elektrophorese	Zelllyse	Zell- Permea- bilisierung	Zellkultur	Pharmazeutische Anwendungen
Brij 35	X		X	X					
CHAPS				X	X				
CTAP			X	X	X		X		
Na-Cholsäure				X					
Na- Deoxycholsäure				X	X				
Digitonin							X		
Na-Docusat									
DM				X					
SDS					X				
L-Sarcosin									
Lubrol						X			
OG					X				
Saponin							X		
SB 10					X				
SB 12				X					
Synperonic F68								X	X
Synperonic F108									X
Triton® X-100		X	X		X				
Tween 20		X		X	X			X	
Tween 80								X	X
ASB-14/16					X				