

Applications (1)	Membrane Proteins/ Proteomic Research	Reconstitution of membrane proteins	Receptor research	Enzymology	Photosynthesis research	Reverse micellar Solubilization	Lipid raft reseach	DNA isolation	RNA isolation
Brij 35	X			X					
CHAPS	X		X	X	X				
CTAP						X		X	X
Na-Cholate		X	X	X					
Na-Deoxycholate			X						
Digitonin	X		X	X	X				
Na-Docusate						X			
DM	X		X	X	X				
SDS	X			X					
L-Sarcosine	X			X				X	X
Lubrol	X			X			X		
OG	X	X	X	X	X				
Saponin									
SB 10	X			X					
SB 12					X				
Synperonic F68									
Synperonic F108									
Triton® X-100	X		X	X	X		X		
Tween 20				X					
Tween 80				X					
ASB-14/16	X								

Applications (2)	Enzyme assays	Immuno assays	Chemiluminescence analysis/Fluorometry	Chromatography	Electrophoresis	Cell lysis	Cell permeabilization	Cell Culture	Pharmaceutical applications
Brij 35	X		X	X					
CHAPS				X	X				
CTAP			X	X	X		X		
Na-Cholate				X					
Na-Deoxycholate				X	X				
Digitonin							X		
Na-Docusate									
DM				X					
SDS					X				
L-Sarcosine									
Lubrol						X			
OG					X				
Saponin							X		
SB 10					X				
SB 12				X					
Synperonic F68								X	X
Synperonic F108									X
Triton® X-100		X	X		X				
Tween 20		X		X	X			X	
Tween 80								X	X
ASB-14/16					X				