

CERTIFICATE OF ANALYSIS

IEF Markers 3-10, SERVA Liquid Mix	Cat.No. : 39212
	Lot No.: P120004

Protein Components	pI	Country of origin (for animals)	Molecular weight (kDa)
Cytochrome C (horse, heart)	10.7	Canada	12.4
Ribonuclease A (bovine pancreas)	9.5	New Zealand /Australia	13.7
Lectin (Lens culinaris)	8.3, 8.0, 7.8		2 x 17.0 / 2 x 8.0
Myoglobin (horse, muscle)	7.4, 6.9	Mexico	17.8
Carbonic anhydrase (bovine erythrocytes)	6.0	USA	29.0
β-Lactoglobulin (bovine, milk)	5.3, 5.2	USA	14.0
Trypsin inhibitor (soybean)	4.5		21.5
Glucose oxidase (Aspergillus niger)	4.2		150.0
Amyloglucosidase (Aspergillus niger)	3.5		97.0
Storage (°C)			+2 to +8
Minimum shelf life			04/ 2013

Handling

ready-to-use

Apply 5µl (50µg protein) per lane when staining with SERVA Violet 17
SERVA Blue R, SERVA Blue G or Coomassie. For silver staining,
dilute 1:10 and apply 5µl (5µg protein) per lane.

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 13.01.2012

Dipl.-Ing. (FH) Bernhard Göckel

Dipl.-Ing. (FH) Christian Nowak

This report has been computer generated and does not contain a signature.

SERVA Electrophoresis GmbH • D-69115 Heidelberg • Carl-Benz-Str. 7 • Tel.: +49 6221 138 40-0 • Fax +49 6221 138 40-10
Geschäftsführer: Dr. Barbara Müller, Dr. Manfred Kurfürst • HRB 336136 AG Mannheim • E-Mail: info.serva.de • http://www.serva.de
COMMERZBANK BLZ: 672 400 39 Kto: 190 38 30 00 • IBAN: DE90 6724 0039 0190 3830 00 • BIC: COBADEFFXXX • Ust.ID-Nr.: DE 812 517 285