

Schneiders Drosophila Pulvermedium, revidiert ohne Calciumchlorid, ohne Natriumbicarbonat	Kat.Nr./Cat.No. :	47521
Schneider's Drosophila powder medium revised without calcium chloride, without sodium bicarbonate		

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		weißes Pulver white powder
Löslichkeit solubility	24.89 g/l in H ₂ O + 0.4 g/l NaHCO ₃ + 0.6 g/l CaCl ₂	klare schwach gelbe Lösung clear slightly yellow solution
Osmolarität (mOsmol/l) osmolarity	24.89 g/l in H ₂ O + 0.4 g/l NaHCO ₃ + 0.6 g/l CaCl ₂ + 1 M NaOH auf/to pH 6.5	270 - 320
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.