

Applications for GEL-FIX™ / NetFix™ Supports

Casting of Gels for Electrophoresis Technique	GEL-FIX™ PAG	GEL-FIX™ Agarose	GEL-FIX™ for Covers	NetFix™ PAG
Native PAGE	yes		yes	yes
IEF, PAGE	yes		yes	yes
IEF, Agarose		yes	yes	
SDS-PAGE	yes		yes	yes
Nucleic acid electrophoresis	yes	yes	yes	yes
Blotting after electrophoresis, PAGE gels			yes	yes
Blotting after electrophoresis, agarose gels			yes	