

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Digitonin wasserlöslich reinst (Digitin) water soluble, research grade	Kat.Nr./Cat.No. : 19551 Lot/Contr.No. : 120308
--	---

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₅₆ H ₉₂ O ₂₉	
Molekulargewicht molecular weight		1 229.34	
Aussehen appearance		weißes Pulver white powder	entspricht corresponds
Wasser (%) water	KF	max. 6.0	1.7
[α]_D/ 20 °C	10 % in 75 % HAc	-49.0° bis/ to -47.0°	-47.9°
Löslichkeit solubility	1 g/ 20 ml H ₂ O	löslich nach Erwärmen auf 95 – 98 °C, bleibt nach Abkühlen für min. 24 Stunden klar dissolved by heating to 95 – 98 °C, remains clear after cooling for at least 24 hours	
Mindesthaltbarkeit minimum shelf life			10/ 2014
Lagerung (°C) storage			+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 20.01.2014

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.