

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Digitonin wasserlöslich reinst (Digitin) water soluble, research grade	Kat.Nr./Cat.No. : 19551 Lot/Contr.No. : 220790
--	---

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₅₆ H ₉₂ O ₂₉	
Molekulargewicht molecular weight		1 229.34	
Aussehen appearance		weißes Pulver white powder	entspricht corresponds
Wasser (%) water	KF	max. 6.0	0.7
[α]_D/ 20 °C	1 % in Ethanol	-57.0° bis/ to -49.0°	-54.7°
Trübung turbidity	1 g wird in 20 ml Wasser suspendiert und durch Erhitzen auf 95 – 98 °C für 15 Minuten gelöst 1 g is suspended in 20 ml of water and dissolved by heating to 95 – 98 °C for 15 minutes		
Mindesthaltbarkeit minimum shelf life			15.11.2024
Lagerung (°C) storage			+15 bis / to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 15.08.2022

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.