

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

| | | |
|---|--------------------------|---------------|
| Glutaraldehyd 25% in Wasser Standardqualität | Kat.Nr./Cat.No. : | 23115 |
| Glutaraldehyde 25% in water standard grade | Lot/Contr.No. : | 240280 |

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|--------------------------|--|---------------------------|
| Summenformel formula | | C ₅ H ₈ O ₂ | |
| Molekulargewicht molecular weight | | 100.13 | |
| Aussehen appearance | | klare, farblose bis gelbliche Lösung clear, colorless to yellowish solution | entspricht corresponds |
| Gehalt (%) assay | | 23.0 – 27.0 | 25.1 |
| Brechungsindex refractive index | 25 °C | 1.3690 – 1.3755 | 1.3734 |
| pH | 25 % in H ₂ O | 3.0 – 4.0 | 3.5 |
| Mindesthaltbarkeit minimum shelf life | | | 18.03.2026 |
| Lagerung (°C) storage | | | +2 bis/ to +8 |

| |
|--|
| <p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p> |
|--|

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 21.03.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.