

CERTIFICATE OF ANALYSIS

Insulin recombinant human	Cat.No.: 26360 Contr.No.: 170986
-------------------------------------	---

Parameter	Method	Specification	Result
Molecular weight		ca. 5800	
Assay (%)	calculated on dry basis	95.0 – 105.0	98.5
Assay (U/mg)	USP calculated on dry basis	min. 27.5	28.5
Assay (U/mg)	as is	N/A	27.1
Loss on drying (%)		max. 10.0	5.0
Bacterial endotoxins (IU/mg)		max. 10	< 0.4
Zinc (%)	calculated on dry basis	0.3 – 0.6	0.39
Minimum shelf life			18.10.2020
Storage (°C)			-15 to -25

The physical/ chemical data corresponds to the requirements of Ph. Eur. and USP.

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 09.11.2017

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.