

**CERTIFICATE OF ANALYSIS**

<b>Myoglobin equine</b> research grade	<b>Cat.No.:</b> : 29895 <b>Contr.No.:</b> : 130087
---	---

Parameter	Method	Specification	Result
<b>Molecular weight</b>		ca. 17 800	
<b>Appearance</b>		red lyophilisate	corresponds
<b>Assay (%)</b>	SDS electrophoresis	min. 95	99
<b>Solubility</b>	1 % in H <sub>2</sub> O	clear red-brown solution	corresponds
<b>Water (%)</b>	KF	max. 10	7
<b>Iron (%)</b>		0.25 – 0.32	0.31
<b>Minimum shelf life</b>			06/ 2015
<b>Storage (°C)</b>			+2 to +8

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date:** 28.01.2013

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Steinmann

This report has been computer generated and does not contain a signature.