

CERTIFICATE OF ANALYSIS

Myoglobin equine research grade	Cat.No.: : 29895 Contr.No.: : 160337
---	---

Parameter	Method	Specification	Result
Molecular weight		ca. 17 800	
Appearance		red lyophilisate	corresponds
Assay (%)	SDS electrophoresis	min. 95	97
Solubility	1 % in H ₂ O	clear red-brown solution	corresponds
Water (%)	KF	max. 10	5
Iron (%)		0.25 – 0.32	0.28
Minimum shelf life			16.11.2019
Storage (°C)			+2 to +8

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 18.04.2016

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

This report has been computer generated and does not contain a signature.