

CERTIFICATE OF ANALYSIS

Myoglobin equine research grade	Cat.No.: : 29895 Contr.No.: : 211025
---	---

Parameter	Method	Specification	Result
Molecular weight		ca. 17 800	
Appearance		red lyophilisate	corresponds
Assay (%)	SDS electrophoresis	min. 95	96
Solubility	1 % in H ₂ O	clear red-brown solution	corresponds
Water (%)	KF	max. 10	9
Iron (%)		0.25 – 0.32	0.30
Minimum shelf life			12.07.2025
Storage (°C)			+2 to +8

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 20.09.2021

Christian Monsler

Marina Konter

This report has been computer generated and does not contain a signature.