

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Nycodenz® 60% (w/v) Lösung in Wasser solution in water	Kat.Nr./Cat.No. : 31000 Lot/Contr.No. : 140845
--	---

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₁₉ H ₂₆ I ₃ N ₃ O ₉	
Molekulargewicht molecular weight		821.1	
Aussehen appearance		klare, farblose Lösung clear, colorless solution	entspricht corresponds
Dichte (g/ml) density	20 °C	1.310 ± 0.002	entspricht corresponds
Brechungsindex refractive index	20 °C	1.4273 ± 0.0003	1.4270
Mindesthaltbarkeit minimum shelf life			08/ 2016
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 21.05.2015

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.