

# CERTIFICATE OF ANALYSIS

<b>Proteinase K from Tritirachium album</b> lyophil.	<b>Cat.No. :</b> 33752 <b>Contr.No.:</b> 240303
---------------------------------------------------------	----------------------------------------------------

Parameter	Method	Specification	Result
<b>Molecular weight</b>		ca. 28 000	
<b>Appearance</b>		white lyophilisate	corresponds
<b>Activity (mAnsonU/mg)</b>		min. 30	45
<b>Extraneous activities</b>	DNase	not detectable	corresponds
	RNase	not detectable	corresponds
<b>Minimum shelf life</b>			01.12.2025
<b>Storage (°C)</b>			+2 to +8

## Unit definition

One unit is defined as the amount of enzyme that liberates Folin-positive amino acids and peptides, corresponding to 1 pmole tyrosine under assay conditions in 1 minute using hemoglobin as substrate.

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date:** 04.04.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.