

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Teepol 610 prakt.	Kat.Nr./Cat.No.:	35796
pract.	Lot/Contr.No. :	121073

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Relative Molmasse relative molecular weight		372	
Aussehen appearance		klare gelbliche bis bräunliche Lösung clear yellowish to brown-colored solution	entspricht corresponds
Anionensidgehalt (%) Anionic active		min. 32.0	32.6
Brechungsindex refractive index	n_D^{20}	1.375 – 1.378	1.376
Identität identity	IR	entspricht corresponds	entspricht corresponds
Formaldehyd (ppm)		1600 – 1950	1877
Dichte (kg/l) density	20°C	1.058 – 1.068	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			02/ 2017
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 26.09.2012

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.