

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>Teepol 610 prakt.</b>	<b>Kat.Nr./Cat.No.:</b>	<b>35796</b>
<b>pract.</b>	<b>Lot/Contr.No. :</b>	<b>131146</b>

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Aussehen</b> appearance		klare gelbliche bis bräunliche Lösung clear yellowish to brown-colored solution	entspricht corresponds
<b>Anionensidagehalt (%)</b> Anionic active		min. 32.0	32.5
<b>Brechungsindex</b> refractive index	$n_D^{20}$	1.375 – 1.378	1.3758
<b>Identität</b> identity	IR	entspricht corresponds	entspricht corresponds
<b>Dichte (kg/l)</b> density	20 °C	1.053 – 1.063	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			10/ 2018
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 17.10.2013**

Dipl.-Ing. (FH) Bernhard Göckel

Jakob Arndt

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.