

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Teepol 610 prakt.	Kat.Nr./Cat.No.:	35796
pract.	Lot/Contr.No. :	140474

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare gelbliche bis bräunliche Lösung clear yellowish to brown-colored solution	entspricht corresponds
Anionensidagehalt (%) Anionic active		min. 32.0	32.5
Brechungsindex refractive index	n_D^{20}	1.375 – 1.378	1.3757
Identität identity	IR	entspricht corresponds	entspricht corresponds
Dichte (kg/l) density	20 °C	1.053 – 1.063	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			04/ 2019
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 15.04.2014

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.