

# CERTIFICATE OF ANALYSIS

**Phosphatase Inhibitor Mix I, powder**

**Cat.No. : 39050**

**Lot No.: P230179**

Components	Target-Phosphatase
Imidazol	Alkaline Phosphatase
Sodium fluoride	Acid Phosphatase Protein-Phosphatase-2A Protein-Phosphatase-2B Protein-Phosphatase-2C
Sodium molybdate	Acid Phosphatase Phosphoprotein Phosphatase
Sodium ortho vanadate	Protein-Tyrosine-Phosphatase Alkaline Phosphatase
Sodium tartrate	Acid Phosphatase

**Storage (°C)** +15 to +30

**Minimum shelf life** 28.04.2025

**Reconstitution** Dissolve the content of the vial (156 mg) in 1 ml distilled water to obtain a 100x stock solution (suitable for treatment of 100 ml cell extract). After reconstitution, store at +2°C to +8°C (stable for 2 years).

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date: 09.05.2023**

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.