

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Natriumacetatpuffer pH 5.2, 3M Lösung	Kat.Nr./Cat.No. :	39572
Sodiumacetate buffer pH 5.2, solution 3M	Lot/Contr.No. :	P140111
molecular biology grade		

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		$C_2H_3O_2Na$	
Molekulargewicht molecular weight		82.0	
Aussehen appearance		klare, farblose Lösung clear, colorless solution	entspricht corresponds
Zusammensetzung composition			
$C_2H_3O_2Na$ (#39571)		246.0 g/l	entspricht corresponds
pH	3 M in H_2O	5.0 – 5.4	5.2
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			02/ 2016
Lagerung (°C) storage			+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 25.02.2014

Dipl.-Ing. (FH) Bernhard Göckel

Jakob Arndt

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.