

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

|                                |                                |
|--------------------------------|--------------------------------|
| <b>TAE Puffer (40x), 1.2 M</b> | <b>Kat.Nr./Cat.No. : 42548</b> |
| TAE buffer (40x)               | <b>Lot/Contr.No. : P160395</b> |

| <b>Parameter</b><br>parameter                   | <b>Spezifikation</b><br>specification               | <b>Ergebnis</b><br>result |
|---|---|---------------------------|
| <b>Aussehen</b><br>appearance                   | klare, farblose Lösung<br>clear, colorless solution | entspricht<br>corresponds |
| <b>pH</b>                                       | 8.3 – 8.7   | 8.4                       |
| <b>Zusammensetzung</b><br>composition           |   |                           |
| TRIS-Acetat                                     | 1.2 M   | entspricht<br>corresponds |
| EDTA-Na <sub>2</sub> ·2H <sub>2</sub> O         | 0.04 M  | entspricht<br>corresponds |
| <b>Mindesthaltbarkeit</b><br>minimum shelf life |   | 15.10.2020                |
| <b>Lagerung (°C)</b><br>storage                 |   | +15 bis/ to +30           |

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
 We do not guarantee that the product can be used for a special application.  
 This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 15.10.2019**

Christian Monsler

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.