

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Sprüh-Magermilchpulver zum Blotten	Kat.Nr./Cat.No. : 42590
Skim milk powder for blotting	Lot/Contr.No. : 230200

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weißes bis gelbliches Pulver white to yellowish powder	entspricht corresponds
Protein (%) protein		32.0 – 40.0	34.1
Fett (%) fat		max. 1,25	0.3
Lactosemonohydrat (%) Lactose monohydrate		45.0 – 56.0	51 – 56
pH	10 % in H ₂ O	6.4 – 6.7	6.7
Asche (%) ash		max. 8.5	7.0
Wassergehalt (%) water content		max. 4.0	3.6
Mindesthaltbarkeit minimum shelf life			12.01.2025
Lagerung (°C) storage			+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 02.03.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.