

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Colcemid Lösung sterilfiltriert	Kat.Nr./Cat.No. :	47253
Colcemid solution sterile filtered	Lot/Contr.No. :	P120033

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, farblose Lösung clear, colorless solution	entspricht corresponds
Konzentration (µg/ml) concentration	in PBS	10	entspricht corresponds
pH		7.0 – 8.0	7.3
Osmolarität (mOsmol/l) osmolarity		max. 300	266
A 1cm/ 246 nm	unverdünnt undiluted	min. 0.8	0.9
A 1cm/ 354 nm	unverdünnt undiluted	min. 0.4	0.4
Sterilität sterility		negativ negative	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			01/ 2014
Lagerung (°C) storage			-15 bis/ to -25

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 30.01.2012

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.