

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Colcemid Lösung sterilfiltriert	Kat.Nr./Cat.No. :	47253
Colcemid solution sterile filtered	Lot/Contr.No. :	P140260

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, farblose Lösung clear, colorless solution	entspricht corresponds
Konzentration (µg/ml) concentration	in PBS	10	entspricht corresponds
pH		7.0 – 8.0	7.4
Osmolarität (mOsmol/l) osmolarity		max. 300	285
A 1cm/ 246 nm	unverdünnt undiluted	min. 0.8	0.955
A 1cm/ 354 nm	unverdünnt undiluted	min. 0.4	0.443
Sterilität sterility		negativ negative	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			06/ 2016
Lagerung (°C) storage			-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 23.06.2014

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.