

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

| | | |
|--|--------------------------|----------------|
| Colcemid Lösung sterilfiltriert | Kat.Nr./Cat.No. : | 47253 |
| Colcemid solution sterile filtered | Lot/Contr.No. : | P150098 |

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|-------------------------|---|---------------------------|
| Aussehen appearance | | klare, farblose Lösung clear, colorless solution | entspricht corresponds |
| Konzentration (µg/ml) concentration | in PBS | 10 | entspricht corresponds |
| pH | | 7.0 – 8.0 | 7.3 |
| Osmolarität (mOsmol/l) osmolarity | | max. 300 | 291 |
| A 1cm/ 246 nm | unverdünnt undiluted | min. 0.8 | 0.961 |
| A 1cm/ 354 nm | unverdünnt undiluted | min. 0.4 | 0.468 |
| Sterilität sterility | | negativ negative | entspricht corresponds |
| Mindesthaltbarkeit minimum shelf life | | | 03/ 2017 |
| Lagerung (°C) storage | | | -15 bis/ to -25 |

| |
|--|
| <p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p> |
|--|

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 16.03.2015

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.