

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Gentamycin Lösung 50mg/ml, sterilfiltriert	Kat.Nr./Cat.No. : 47991
Gentamycin solution 50mg/ml, sterile filtered	Lot/Contr.No. : P160013

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, leicht gelbliche Lösung clear, slightly yellowish solution	entspricht corresponds
Konzentration (mg/ ml) concentration	in H ₂ O	80 mg Gentamycinsulfat entspricht 50 mg Gentamycin 80 mg Gentamycin sulfate corresponds 50 mg Gentamycin	
pH		4.0 – 6.0	4.7
Osmolarität (mOsmol/ l) osmolarity		ca. 100	132
Sterilität sterility		negativ negative	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			14.01.2018
Lagerung (°C) storage			-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH **Druckdatum/ printing date: 25.01.2016**
Qualitätskontrolle/ Quality Control

Dipl.-Ing. (FH) Bernhard Göckel Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.