

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Gentamycin Lösung 50mg/ml, sterilfiltriert	Kat.Nr./Cat.No. : 47991
Gentamycin solution 50mg/ml, sterile filtered	Lot/Contr.No. : P210109

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, leicht gelbliche Lösung clear, slightly yellowish solution	entspricht corresponds
Konzentration (mg/ ml) concentration	in H ₂ O	80 mg Gentamycinsulfat entspricht 50 mg Gentamycin 80 mg Gentamycin sulfate corresponds 50 mg Gentamycin	
pH		4.0 – 6.0	5.0
Osmolarität (mOsmol/ l) osmolarity		ca. 100	141
Sterilität sterility		negativ negative	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			06.04.2025
Lagerung (°C) storage			-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 14.03.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.