

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Gentamycin Lösung 50mg/ml, sterilfiltriert	Kat.Nr./Cat.No. : 47991
Gentamycin solution 50mg/ml, sterile filtered	Lot/Contr.No. : P230379

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, leicht gelbliche Lösung clear, slightly yellowish solution	entspricht corresponds
Konzentration (mg/ ml) concentration	in H ₂ O	80 mg Gentamycinsulfat entspricht 50 mg Gentamycin 80 mg Gentamycin sulfate corresponds 50 mg Gentamycin	
pH		4.0 – 6.0	4.6
Osmolarität (mOsmol/l) osmolarity		ca. 100	120
Sterilität sterility		negativ negative	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			23.10.2025
Lagerung (°C) storage			-15 bis/ to -25

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH **Druckdatum/ printing date: 31.10.2023**
Qualitätskontrolle/ Quality Control

Dr. Judith Koch Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.