

**ANALYSENZERTIFIKAT**  
 CERTIFICATE OF ANALYSIS

<b>Acrylamid 2x, reinst</b>	<b>Kat.Nr./Cat.No. :</b>	<b>10675</b>
Acrylamide 2x, research grade	<b>Lot/Contr.No. :</b>	<b>190315</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Summenformel</b> formula		C <sub>3</sub> H <sub>5</sub> NO	
<b>Molekulargewicht</b> molecular weight		71.1	
<b>Aussehen</b> appearance		weiße Kristalle white crystals	entspricht corresponds
<b>Reinheit (%)</b> purity	Titration	min. 98.0	99.96
<b>Acrylsäure (%)</b> acrylic acid	Titration	max. 0.03	< 0.01
<b>A 1 cm/ 290 nm</b>	5% in H <sub>2</sub> O	max. 0.7	0.442
<b>Löslichkeit</b> solubility	5% in H <sub>2</sub> O	klar und farblos clear and colorless	entspricht corresponds
<b>pH</b>	5% in H <sub>2</sub> O	5.0 – 8.0	6.1
<b>Leitfähigkeit (µS/cm)</b> conductivity	40% in H <sub>2</sub> O	max. 100	< 1
<b>Polymerisationszeit (min)</b> polymerizing time	15% Gel, 3 mM APS/ TEMED/ 25 °C	max. 30	15
<b>Mindesthaltbarkeit</b> minimum shelf life			25.01.2024
<b>Lagerung (°C)</b> storage			+2 bis/ to +8

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b>  <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b>          We do not guarantee that the product can be used for a special application.          This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**
**Druckdatum/ printing date:** 26.03.2019

Christian Monsler

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.