

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Agar Agar SERVA hohe Gelstärke Pulver reinst **Kat.Nr./Cat.No. : 11396**
Agar Agar SERVA high gel-strength powder research grade **Lot/Contr.No. : 120264**

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		hellelfenbeinfarbenes Pulver light ivory-colored powder	entspricht corresponds
Gelstärke (g/cm²) gel strength	1.5 % Gel 1.5 % gel	min. 700	1000
Gelertemperatur (°C) point of solidification	1.5 % Lösung 1.5 % solution	32 – 35	35
A 1cm/430nm	1.5 % in H ₂ O (60 °C)	max. 0.15	0.098
A 1cm/525nm	1.5 % in H ₂ O (60 °C)	max. 0.06	0.041
Trocknungsverlust (%) loss on drying	4 h, 105 °C	max. 10.0	8.1
Asche (%) ash		max. 4.0	3.4
pH	1.5 % in H ₂ O (60 °C)	6.0 – 8.0	7.0
Mindesthaltbarkeit minimum shelf life			03/ 2017
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 20.03.2012

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.