

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Agar Agar SERVA hohe Gelstärke Pulver reinst	Kat.Nr./Cat.No. :	11396
Agar Agar SERVA high gel-strength powder research grade	Lot/Contr.No. :	160661

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		hell elfenbeinfarbenes bis beiges Pulver light ivory-colored to beige powder	entspricht corresponds
Gelstärke (g/cm²) gel strength	1.5 % Gel 1.5 % gel	min. 700	900
Gelertemperatur (°C) point of solidification	1.5 % Lösung 1.5 % solution	32 – 35	35
A 1cm/430nm	1.5 % in H ₂ O (60 °C)	max. 0.2	0.190
A 1cm/525nm	1.5 % in H ₂ O (60 °C)	max. 0.1	0.099
Trocknungsverlust (%) loss on drying	4 h, 105 °C	max. 10.0	6.8
Asche (%) ash		max. 4.0	3.4
pH	1.5 % in H ₂ O (60 °C)	6.0 – 8.0	7.4
Mindesthaltbarkeit minimum shelf life			23.06.2020
Lagerung (°C) storage			+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 06.07.2016

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.