

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

<b>Agarose SERVA neutral für IEF</b>	<b>Kat.Nr./Cat.No. : 11402</b>
Agarose SERVA neutral for IEF	<b>Lot/Contr.No. : 190390</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		weißes Pulver white powder	entspricht corresponds
<b>Gelstärke (g/cm<sup>2</sup>)</b> gel strength	1.5 % Gel 1.5 % gel	> 800	1 200
<b>Trocknungsverlust (%)</b> loss on drying	4h, 105°C	max. 5.0	3.3
<b>Elektroendosmose (mr)</b> electroendosmosis		± 0.0	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			17.04.2024
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b>  <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b>          We do not guarantee that the product can be used for a special application.          This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 26.04.2019**

Christian Monsler

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.