

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Albumin aus Hühnerei, lyophil.	Kat.Nr./Cat.No. : 11842
Albumin egg (ovalbumin) lyophil.	Lot/Contr.No. : 180080

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weißes bis hellgelbes Lyophilisat white to pale yellow lyophilisate	entspricht corresponds
Reinheit (Ovalbumin) (%) purity (ovalbumin)	SDS-PAGE	min. 90	> 90
Gehalt (Protein) (%) assay (protein)		min. 85	93
Asche (%) ash		max. 5.0	< 0.1
Trocknungsverlust (%) loss on drying		max 6.0	1.8
Ursprung origin		Hühnereiweiß chicken egg white	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			17.01.2021
Lagerung (°C) storage			-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 22.10.2019

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.