

CERTIFICATE OF ANALYSIS

Albumin human lyophilized	Cat.No. : 11877 Contr.No.: 130189
-------------------------------------	--

Parameter	Method	Specification	Result
Appearance		pale yellow to greenish flakes	corresponds
Albumin purity (%)		> 97	99
Total Proteins (%)		> 92	95
Moisture (%)		< 5	1.8
pH	10 % in H ₂ O	6.5 – 7.4	6.8
Mesophilic germs (CFU/ g)		< 100	corresponds
Minimum shelf life			18.07.2018
Storage (°C)			+2 to +8

The source material was tested and found negative for:

- Antibodies: anti-HIV 1 + 2, HCV, Hbc, HTLV 1 + 2
- HBs antigen
- HIV 1 and HCV RNA
- Syphilitic serology

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 26.04.2017

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.