

CERTIFICATE OF ANALYSIS

Albumin human lyophilized	Cat.No. : 11877 Contr.No.: 150627
-------------------------------------	--

Parameter	Method	Specification	Result
Appearance		pale yellow to greenish flakes	corresponds
Albumin purity (%)		> 97	97
Total Proteins (%)		> 92	96
Moisture (%)		< 5	2.1
pH	10 % in H ₂ O	6.5 – 7.4	6.9
Mesophilic germs (CFU/ g)		< 100	corresponds
Minimum shelf life			01/ 2018
Storage (°C)			+2 to +8

The source material was tested and found negative for:

- Antibodies: anti-HIV 1 + 2, anti-HCV
- HBs antigen
- HIV 1 and HCV RNA
- Syphilitic serology

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 02.06.2015

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

This report has been computer generated and does not contain a signature.