

CERTIFICATE OF ANALYSIS

Albumin bovine Fraction V	Cat.No. : 11945
Protease- and Fatty-Acid-free, diagnostic grade, lyophil.	Contr.No.: 151372

Parameter	Method	Specification	Result
Molecular weight		ca. 67 000	
Appearance	Visual inspection	white to yellow with tan to green cast lyophilisate	corresponds
Solubility	4 % in H ₂ O	clear to slightly hazy, faint yellow to yellow-green solution	corresponds
Protein (%) protein factor 6.22	Nitrogen analyzer dry weight basis	≥ 98	101.9
Purity (%)	Agarose zone electrophoresis	≥ 98	100
pH	10 % solution	6.5 – 7.5	6.9
Moisture (%)	Loss on drying	0 – 5.0	0.9
Ash (%)	Residue on ignition	≤ 2	0.4
IgG (bovine)	Turbidimetric analysis	not detected	not detected
Protease	Enzymatic assay	≤ 0.005 U/mg	< 0.005 U/mg
Fatty acids (%)	Enzymatic assay	< 0.01	< 0.001
Heavy metals (ppm)	Atomic absorption	< 10	< 0.50
Endotoxin (EU/ mg)	LAL	< 3	< 0.1
Minimum shelf life			09.04.2019
Storage (°C)			+2 to +8

This product was derived from bovine blood collected at a USDA licensed establishment located in the United States. The animals received ante- and post-mortem health inspection at the abattoir by a Veterinary Medical Officer, USDA Food Safety and Inspection Services, and they were apparently free from infectious and contagious diseases and injuries parasites.

During the manufacturing process, the material was subjected to a pH of ≤5 and a temperature in excess of 65°C for at least 3 hours and was not commingled with any other material of animal origin. Blood collection records are incorporated into the manufacturing batch documents.

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Quality Control

Printing date: 04.11.2015

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.