

CERTIFICATE OF ANALYSIS

Albumin bovine Fraction V, pH 7.0 Life Science grade, lyophil.	Cat.No. : 11946 Contr.No.: 091228
--	--

Parameter	Method	Specification	Result
Molecular Weight		ca. 67 000	
Appearance	Visual Inspection	white to yellow with tan to green cast lyophilisate	corresponds
Protein (%) (protein factor 6.25)	Nitrogen Analyzer, Anhydrous	≥ 96	103
Purity (%) (Albumin, 7 % solution)	Cellulose Acetate Electrophoresis	≥ 98	100
pH	1 % solution in 0.15 M NaCl	6.8 – 7.2	7.0
Moisture (%)	Loss on Drying	≤ 5.0	1.9
Sulphated Ash (%)	Residue on Ignition	≤ 2.0	1.2
IgG (µg/g)	By ELISA	≤ 50	1.0
Solubility	4 % Solution @ A550	≤ 0.1	< 0.1
Color	4 % Solution @ A404	≤ 0.15	0.01
Heavy Metals (ppm)	By USP	≤ 20	< 10
Mycoplasma		None detected	None detected
Minimum shelf life			05/ 2014
Storage (°C)			+2 to +8

The product has been subjected to a heat treatment at a room temperature of at least 65 °C for a period of at least 3 hours and subjected to pH 5 or lower for a period of at least 3 hours.

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>

SERVA Electrophoresis GmbH
Quality Control

Printing date: 10.04.2012

Dipl.-Ing. (FH) Bernhard Göckel

Dipl.-Ing. (FH) Christian Nowak

This report has been computer generated and does not contain a signature.