

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Bathocuproin disulfonsäure·Na₂-salz	p.A.	Kat.Nr./Cat.No. :	14470
Bathocuproine disulfonic acid·Na ₂ -salz	analytical grade	Lot/Contr.No. :	180513

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₂₆ H ₁₈ N ₂ O ₆ S ₂ ·Na ₂	
Molekulargewicht molecular weight		564.6	
Aussehen appearance		beiges Pulver beige powder	entspricht corresponds
Gehalt (%) assay	Titration	min. 98.0	100.8
Wasser (%) water	KF	max. 4.0	1.5
Identität identity	IR	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			07.05.2025
Lagerung (°C) storage			+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control
Druckdatum/ printing date: 19.04.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.