

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Bathocuproin disulfonsäure·Na₂-salz	p.A.	Kat.Nr./Cat.No. :	14470
Bathocuproine disulfonic acid·Na ₂ -salz	analytical grade	Lot/Contr.No. :	221090

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₂₆ H ₁₈ N ₂ O ₆ S ₂ ·Na ₂	
Molekulargewicht molecular weight		564.6	
Aussehen appearance		beiges Pulver beige powder	entspricht corresponds
Gehalt (%) assay	Titration	min. 98.0	100.7
Wasser (%) water	KF	max. 4.0	3.6
Identität identity	IR	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			09.11.2027
Lagerung (°C) storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 15.11.2022

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.