

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Brij 35 pract.	Kat.Nr./Cat.No. :	15230
Polyoxyethylene monolauryl ether; n ca. 23	Lot/Contr.No. :	180640

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Molekulargewicht molecular weight		ca. 1 200	
Aussehen appearance		weiße, wachsartige Schuppen white, wax-like flakes	entspricht corresponds
Wasser (%) water		max. 1.0	1
Hydroxylzahl (mg KOH/ g) hydroxyl value		40 – 60	43
Säurezahl (mg KOH/ g) acid value		max. 1.0	0.1
Verseifungszahl (mg KOH/ g) saponification value		max. 3.0	0
Mindesthaltbarkeit minimum shelf life			19.05.2020
Lagerung (°C) storage			+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 18.06.2018

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.