

CERTIFICATE OF ANALYSIS

Carbonic anhydrase from bovine erythrocytes EC 4.2.1.1 lyophil. salt-free	Cat.No. : 15882 Contr.No.: 150724
---	--

Parameter	Method	Specification	Result
Molecular weight		ca. 29 000	
Appearance		white lyophilisate	corresponds
Solubility	1 mg/ml 0.04 M Tris/HCl- buffer, pH 7.6	clear and colorless	corresponds
Activity (U/mg)	Carbonic anhydrase	min. 1.0	1.1
Protein concentration (g protein/ g material)	Lowry		0.94
Minimum shelf life			02.05.2018
Storage (°C)			-15 to -25

Unit definition

1 Unit is that amount of enzymatic activity that catalyzes the hydrolysis of 1 µmol of p-nitrophenylacetate per minute at 25°C, pH 7.6.

We hereby certify that the raw material used was of Australian origin, collected from abattoirs which are under veterinary control and from healthy animals certified as fit for human consumption. All reasonable precautions were taken to prevent contamination of the products during and after processing. The above product is supplied for research purposes and/or in vitro use only, and is not to be used in humans or animals.

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>

SERVA Electrophoresis GmbH
Quality Control

Printing date: 26.04.2017

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.