

CERTIFICATE OF ANALYSIS

Carbonic anhydrase from bovine erythrocytes	Cat.No. : 15882
EC 4.2.1.1 lyophil. salt-free	Contr.No.: 190886

Parameter	Method	Specification	Result
Molecular weight		ca. 29 000	
Appearance		white lyophilisate	corresponds
Solubility	1 mg/ml 0.04 M Tris/HCl- buffer, pH 7.6	clear and colorless	corresponds
Activity (U/mg)	Carbonic anhydrase	min. 1.0	1.2
Protein concentration (g protein/ g material)	Lowry		0.830
Minimum shelf life			05.09.2021
Storage (°C)			-15 to -25

Unit definition

1 Unit is that amount of enzymatic activity that catalyzes the hydrolysis of 1 µmol of p-nitrophenylacetate per minute at 25°C, pH 7.6.

We hereby certify that the raw material used was of Australian origin, collected from abattoirs which are under veterinary control and from healthy animals certified as fit for human consumption. All reasonable precautions were taken to prevent contamination of the products during and after processing. The above product is supplied for research purposes and/or in vitro use only, and is not to be used in humans or animals.

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control
Printing date: 05.09.2019

Christian Monsler

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.