

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Dextran 4	Kat.Nr./Cat.No. :	18687
technical grade	Lot/Contr.No. :	230317

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weißes bis schwach grauweißes Pulver white to slightly off-white powder	entspricht corresponds
Trocknungsverlust (%) loss on drying		max. 7.0	4.4
pH	10 % in H ₂ O	4.5 – 7.0	4.8
Mindesthaltbarkeit minimum shelf life			13.03.2028
Lagerung (°C) storage			+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 20.03.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.