

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Dextran 500 techn. grade	Kat.Nr./Cat.No. :	18696
	Lot/Contr.No. :	120148

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weißes bis grauweißes Pulver white to off-white powder	entspricht corresponds
Mittleres Molekulargewicht average molecular weight	Gelfiltration	ca. 500 000	entspricht corresponds
Trocknungsverlust (%) loss on drying		max. 7.0	5.1
pH	6 % in H ₂ O	4.5 – 7.0	6.2
A 1cm/ 375 nm	6 % in H ₂ O	max. 0.15	0.05
Mindesthaltbarkeit minimum shelf life			11/ 2015
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 09.02.2012

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.