

ANALYSENZERTIFIKAT
 CERTIFICATE OF ANALYSIS

Dextran 500 techn. grade	Kat.Nr./Cat.No. : 18696
	Lot/Contr.No. : 140386

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weißes bis grauweißes Pulver white to off-white powder	entspricht corresponds
Mittleres Molekulargewicht average molecular weight	Gelfiltration	ca. 500 000	entspricht corresponds
Trocknungsverlust (%) loss on drying		max. 7.0	5.7
pH	6 % in H ₂ O	4.5 – 7.0	5.5
A 1cm/ 375 nm	6 % in H ₂ O	max. 0.15	0.02
Mindesthaltbarkeit minimum shelf life			30.11.2019
Lagerung (°C) storage			+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 13.06.2018

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.