

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Digitonin wasserlöslich reinst (Digitin) water soluble, research grade	Kat.Nr./Cat.No. : 19551 Lot/Contr.No. : 150896
--	---

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₅₆ H ₉₂ O ₂₉	
Molekulargewicht molecular weight		1 229.34	
Aussehen appearance		weißes Pulver white powder	entspricht corresponds
Wasser (%) water	KF	max. 6.0	0.9
[α]_D/ 20 °C	10 % in 75 % HAc	-50.0° bis/ to -47.0°	-48.9°
Trübung turbidity	1 g wird in 20 ml Wasser suspendiert und durch Erhitzen auf 95 – 98 °C für 15 Minuten gelöst 1 g is suspended in 20 ml of water and dissolved by heating to 95 – 98 °C for 15 minutes		
Mindesthaltbarkeit minimum shelf life			02.12.2017
Lagerung (°C) storage			+15 bis / to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 07.03.2017

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.