

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

2-Dodecenylnbernsteinsäureanhydrid pract.	Kat.Nr./Cat.No. :	20755
2-Dodecenylnsuccinic acid anhydride pract.	Lot/Contr.No. :	150714

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₁₆ H ₂₆ O ₃	
Molekulargewicht molecular weight		266.38	
Aussehen appearance		klare, leicht gelbe, viskose Flüssigkeit clear, light yellow, viscous liquid	entspricht corresponds
Identität identity	IR	entspricht corresponds	entspricht corresponds
Farbzahl color	APHA	max. 400	400
Brechungsindex refractive index	n _D ^{20°C}	1.470 – 1.478	1.478
Mindesthaltbarkeit minimum shelf life			03.03.2019
Lagerung (°C) storage			+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 07.12.2017

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.