

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

2-Dodecenylnbernsteinsäureanhydrid pract.	Kat.Nr./Cat.No. :	20755
2-Dodecenylnsuccinic acid anhydride pract.	Lot/Contr.No. :	230950

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₁₆ H ₂₆ O ₃	
Molekulargewicht molecular weight		266.38	
Aussehen appearance		klare, leicht gelbe, viskose Flüssigkeit clear, light yellow, viscous liquid	entspricht corresponds
Identität identity	IR	entspricht corresponds	entspricht corresponds
Farbzahl color	Gardner	max. 6	5
Brechungsindex refractive index	n _D ^{20°C}	1.470 – 1.478	1.478
Mindesthaltbarkeit minimum shelf life			14.09.2025
Lagerung (°C) storage			+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 12.10.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.